
shampoo Documentation

Release 0.1.dev156

Brett Morris

November 30, 2016

I	Links	3
II	General Documentation	7
1	Installation	9
1.1	Requirements	9
1.2	Install shampoo	9
1.3	Testing	10
2	Getting Started	11
2.1	Simple numerical reconstruction	11
3	Tutorials	13
3.1	Reconstruction Tutorials	13
4	Reference/API	17
4.1	shampoo Package	17
III	Authors	29
	Bibliography	33
	Python Module Index	35

shampoo is an open source Python package for digital holographic microscopy with [SHAMU](#), the Submersible Holographic Astrobiology Microscope with Ultraresolution. At present, shampoo does holographic reconstruction, and in the near future, shampoo will support specimen detection and tracking.

Warning: This code is in development and may break without warning.

Part I

Links

- [Source code](#)
- [Issues](#)
- [Docs](#)

Part II

General Documentation

Installation

1.1 Requirements

Note: Users are strongly recommended to manage these dependencies with the excellent [Anaconda Python Distribution](#) which provides easy access to all of the above dependencies and more.

shampoo works on Linux, Mac OS X and Windows. It requires Python 3.5 or 2.7 (earlier versions are not supported) as well as the following packages:

- [Numpy](#)
- [scipy](#)
- [Matplotlib](#)
- [skimage](#)
- [sklearn](#)
- [Astropy](#)
- [h5py](#)
- [pyfftw](#)

1.1.1 pyFFTW

shampoo depends on a package called [pyfftw](#) for speedy, multithreaded Fourier transforms, which is easy to install on Mac OS X and linux but may be tricky on Windows machines. We recommend that Windows users install pyfftw by doing the following steps via conda:

```
conda install -c salilab fftw  
pip install pyfftw
```

1.2 Install shampoo

You can install the latest developer version of shampoo by cloning the git repository:

```
git clone https://github.com/bmorris3/shampoo.git
```

...then installing the package with:

```
cd shampoo
python setup.py install
```

1.3 Testing

If you want to check that all the tests are running correctly with your Python configuration, start up python, and type:

```
import shampoo
shampoo.test()
```

If there are no errors, you are good to go!

1.3.1 More

shampoo follows [Astropy](#)'s guidelines for affiliated packages—installation and testing for the two are quite similar! Please see [Astropy's installation page](#) for more information.

Getting Started

2.1 Simple numerical reconstruction

shampoo comes with a sample hologram of a US Air Force resolution target to reconstruct in the data directory. Let's reconstruct it with shampoo.

Let's start in the shampoo repository's top level directory. First, one must import shampoo, and specify the path to the file to reconstruct and the propagation distance:

```
from shampoo import Hologram
hologram_path = 'data/USAF_test.tif'
propagation_distance = 0.03685 # m
```

Now we will import the hologram, which is stored in TIF format, using the `from_tif` method:

```
h = Hologram.from_tif(hologram_path)
```

If your digital holographic microscope has different properties from SHAMU, for example, a different wavelength and pixel size, you will want to set those properties in the `Hologram` constructor, like this:

```
h = Hologram.from_tif(hologram_path, wavelength=650e-6, dx=3.5e-6, dy=3.5e-6)
```

The defaults are set for SHAMU's configuration, so the defaults will work for this example hologram which was recorded by SHAMU.


Now let's reconstruct the hologram:

```
wave = h.reconstruct(propagation_distance)
```

That's it! The reconstructed wave is stored within `wave`, which is a `ReconstructedWave` object that stores the 2D complex `ndarray` reconstructed wave, plus has attributes to get the phase and intensity arrays (`phase` and `intensity`), and a plotting function `plot` which we will use to plot the phase and intensity from the reconstructed wave:

```
import matplotlib.pyplot as plt
fig, ax = wave.plot()
fig.suptitle("USAF Target")
fig.tight_layout()
plt.show()
```

Here's the result:



Now you're ready to reconstruct your holograms!

Warning: This release should be considered a “preview”, as shampoo is still under development.

For more tutorials, see the Tutorials documentation.

[Return to Top](#)

Tutorials

We've put together some basic tutorials to show you how shampoo works.

If you're looking for details about specific functions' inputs/outputs, see instead the [Reference/API](#).

If you have any trouble reproducing these examples, you are welcome to submit an issue on our GitHub [issue tracker](#), or to submit a [pull request](#) with a fix, if you're adventurous.

We currently have the following tutorials:

3.1 Reconstruction Tutorials

If you're looking for a quick reference for beginning your first reconstruction, see [Getting Started](#). Below we'll dive into some more in-depth examples.

3.1.1 Contents

- *Reconstructing multiple z-planes*
- *Reconstructing a cropped hologram*

3.1.2 Reconstructing multiple z-planes

If you want to reconstruct several z-planes – for example, to find the propagation distance to the USAF test target – you can use the following pattern. First, we will load the raw hologram (only once):

```
from shampoo import Hologram
import numpy as np

hologram_path = 'data/USAF_test.tif'
h = Hologram.from_tif(hologram_path)
```

Next we must set the range of propagation distances to reconstruct. We will use `linspace` to create a linearly-spaced array of five propagation distances:

```
n_z_slices = 5
propagation_distances = np.linspace(0.03585, 0.03785, n_z_slices)
```

Then we will loop over the propagation distances, calling `reconstruct` for each one. We will store the reconstructed wave intensities into a data cube called `intensities`:

```
# Allocate some memory for the complex reconstructed waves
intensities = np.zeros((n_z_slices, h.hologram.shape[0], h.hologram.shape[1]),
                      dtype=np.complex128)

# Loop over all propagation distances
for i, distance in enumerate(propagation_distances):

    # Reconstruct at each distance
    wave = h.reconstruct(distance)


    # Save the result to the data cube
    intensities[i, ...] = wave.intensity
```

Now `intensities` contains all of the reconstructed wave intensities, so if we want to find the propagation distance at best focus, we could apply a very crude focus metric to the intensity arrays. For example, the standard deviation of the intensities will be maximal at the propagation distance nearest to the USAF target, so let's plot the standard deviation of the intensities at each propagation distance:

```
# Measure standard deviation within each intensity image
std_intensities = np.std(intensities, axis=(1, 2))

# Initialize a figure object
import matplotlib.pyplot as plt

fig, ax = plt.subplots()
ax.plot(propagation_distances, std_intensities, 'o-')
ax.axvline(propagation_distances[np.argmax(std_intensities)], ls='--',
           label='Best propagation distance')
ax.set(xlabel='Propagation distance', ylabel='std( Intensity )')
ax.legend()
plt.show()
```



You can see that the best propagation distance is near the same distance that we used in the first example (that's why we picked it!).

[Return to Top](#)

3.1.3 Reconstructing a cropped hologram

Sometimes you need speed. Sometimes you can afford to reconstruct only the central part of the hologram's field of view if it will save you some time. `Hologram` has a built-in option to help you in this situation.

The keyword argument `crop_fraction` sets the fraction of the original image to crop out. By default it is set to `None`. If you set `crop_fraction=0.5`, a 1024x1024 hologram will be cropped to 512x512 before being reconstructed:


```
# Import package, set hologram path, propagation distance
from shampoo import Hologram
hologram_path = 'data/USAF_test.tif'
propagation_distance = 0.03685 # m

# Construct the hologram object, reconstruct the complex wave
h = Hologram.from_tif(hologram_path, crop_fraction=0.5)
wave = h.reconstruct(propagation_distance)

# Plot the reconstructed phase/intensity
import matplotlib.pyplot as plt
fig, ax = wave.plot()
```

```
fig.suptitle("USAF Target")  
fig.tight_layout()  
plt.show()
```

Since the implementations of the FFT used by shampoo generally run on order $O(N \log N)$ time, the cropped hologram reconstruction should compute about 2.2x faster than the full-sized one.



[Return to Top](#)

Reference/API

4.1 shampoo Package

4.1.1 Functions

<code>cluster_focus_peaks(xyz[, eps, min_samples])</code>	Use DBSCAN to identify single particles through multiple focus planes.
<code>create_hdf5_archive(hdf5_path, ...[, ...])</code>	Create HDF5 file structure for holograms and phase/intensity reconstructions.
<code>fftshift(x[, additional_shift, axes])</code>	Shift the zero-frequency component to the center of the spectrum, or with some
<code>find_focus_plane(roi_cube[, focus_on, plot])</code>	Find focus plane in a cube of reconstructed waves at different propagation dista
<code>glue_focus(xyz, labels)</code>	Launch a glue session with focusing results.
<code>locate_specimens(wave_cube, positions, ...)</code>	Identify the (x, y, z) coordinates of a specimen.
<code>open_hdf5_archive(hdf5_path)</code>	Load and return a shampoo HDF5 archive.
<code>save_scaled_image(image, filename[, margin, ...])</code>	Save an image to png.
<code>test([package, test_path, args, plugins, ...])</code>	Run the tests using <code>py.test</code> .
<code>unwrap_phase(reconstructed_wave[, seed])</code>	2D phase unwrap a complex reconstructed wave.

`cluster_focus_peaks`

`shampoo.cluster_focus_peaks(xyz, eps=5, min_samples=3)`

Use DBSCAN to identify single particles through multiple focus planes.

Parameters

xyz : `ndarray`

Matrix of (x, y, z) positions for each peak detected

eps : float

Passed to the `eps` argument of DBSCAN

min_samples : int

Passed to the `min_samples` argument of DBSCAN

Returns

labels : `ndarray`

List of cluster labels for each peak. Labels of -1 signify noise points.

create_hdf5_archive

`shampoo.create_hdf5_archive(hdf5_path, hologram_paths, n_z, metadata={}, compression='lzf', overwrite=False)`

Create HDF5 file structure for holograms and phase/intensity reconstructions.

Parameters

hdf5_path : string

Name of new HDF5 archive

hologram_paths : string

List of all holograms

n_z : int

Number of z-stacks to allocate space for

meta : dict

Metadata to store with in top-level of the HDF5 archive

Returns

f : File

Opened HDF5 file

fftshift

`shampoo.fftshift(x, additional_shift=None, axes=None)`

Shift the zero-frequency component to the center of the spectrum, or with some additional offset from the center.

This is a more generic fork of `fftshift`, which doesn't support additional shifts.

Parameters

x : array_like

Input array.

additional_shift : list of length M

Desired additional shifts in x and y directions respectively

axes : int or shape tuple, optional

Axes over which to shift. Default is None, which shifts all axes.

Returns

y : ndarray

The shifted array.

find_focus_plane

`shampoo.find_focus_plane(roi_cube, focus_on='amplitude', plot=False)`

Find focus plane in a cube of reconstructed waves at different propagation distances.

Uses the autofocus method of Dubois et al. 2006 who showed that the integral over the image plane of the amplitude of the reconstructed wave is minimum at the focal plane for a pure amplitude object [R4]. This will only work for small cubes centered on the specimen.

Parameters**roi_cube** : `ndarray`

Reconstructed waves at N propagation distances with M by M pixels, with a shape of (N, M, M)

focus_on : {"amplitude", "phase"} (optional)

Focus on the phase or amplitude?

plot : bool (optional)

Make plots of the integral of the amplitude of the reconstructed wave as a function of distance. Default is False.

Returns**focus_index** : int

Index of the z-plane that is in focus

significance : float

Significance of a detection of a particle in the ROI cube, measured by the significance of the negative peak in the derivative of the sum of the unwrapped phase in the image plane, with respect to propagation distance. For example, if `significance < 3` the detection of a specimen may be legitimate.

glue_focus`shampoo.glue_focus(xyz, labels)`

Launch a glue session with focusing results.

Parameters**xyz** : `ndarray`

Matrix of (x, y, z) positions of particles

labels : `ndarray`

Labels for particles assigned by clustering algorithm

Returns**ga** : `GlueApplication`

Glue qt GUI application session

locate_specimens`shampoo.locate_specimens(wave_cube, positions, labels, distances, plots=False)`

Identify the (x, y, z) coordinates of a specimen.

Parameters**wave_cube** : `ndarray`

Cube of reconstructed waves

positions : `ndarray`

(x,y,z) positions of objects detected by the blob finder

labels : `ndarray`

Clustering labels for each (x,y,z) coordinate, identifying groups of positions, i.e., single particles detected at multiple z-planes

distances : `ndarray`

Propagation distances, same length as the first axis of `complex_cube`

Returns

specimen_coordinates : `ndarray`

(x, y, z) coordinates of each detected specimen

specimen_coordinates : `ndarray`

Significance of each specimen detection. See docs of `find_focus_plane` for hints on how to interpret the significance quantity.

open_hdf5_archive

`shampoo.open_hdf5_archive(hdf5_path)`

Load and return a shampoo HDF5 archive.

Parameters

hdf5_path : string

Name of HDF5 archive

Returns

f : File

Opened HDF5 file

save_scaled_image

`shampoo.save_scaled_image(image, filename, margin=100, blobs=None, min=0.01, max=99.99)`

Save an image to png.

Parameters

image : `ndarray`

Image to save

filename : str

Path to where to save the png file

blobs : list or `ndarray` or `None`

(x, y, z) positions

margin :

min : float

Colormap scaling minimum

max : float

Colormap scaling maximum

test

shampoo.**test**(*package=None, test_path=None, args=None, plugins=None, verbose=False, pastebin=None, remote_data=False, pep8=False, pdb=False, coverage=False, open_files=False, **kwargs*)
 Run the tests using `py.test`. A proper set of arguments is constructed and passed to `pytest.main`.

Parameters

package : str, optional

The name of a specific package to test, e.g. 'io.fits' or 'utils'. If nothing is specified all default tests are run.

test_path : str, optional

Specify location to test by path. May be a single file or directory. Must be specified absolutely or relative to the calling directory.

args : str, optional

Additional arguments to be passed to `pytest.main` in the `args` keyword argument.

plugins : list, optional

Plugins to be passed to `pytest.main` in the `plugins` keyword argument.

verbose : bool, optional

Convenience option to turn on verbose output from `py.test`. Passing True is the same as specifying '-v' in args.

pastebin : {'failed', 'all', None}, optional

Convenience option for turning on `py.test` pastebin output. Set to 'failed' to upload info for failed tests, or 'all' to upload info for all tests.

remote_data : bool, optional

Controls whether to run tests marked with `@remote_data`. These tests use online data and are not run by default. Set to True to run these tests.

pep8 : bool, optional

Turn on PEP8 checking via the `pytest-pep8` plugin and disable normal tests. Same as specifying '--pep8 -k pep8' in args.

pdb : bool, optional

Turn on PDB post-mortem analysis for failing tests. Same as specifying '--pdb' in args.

coverage : bool, optional

Generate a test coverage report. The result will be placed in the directory `htmlcov`.

open_files : bool, optional

Fail when any tests leave files open. Off by default, because this adds extra run time to the test suite. Requires the `psutil` package.

parallel : int, optional

When provided, run the tests in parallel on the specified number of CPUs. If parallel is negative, it will use the all the cores on the machine. Requires the `pytest-xdist` plugin installed. Only available when using Astropy 0.3 or later.

kwargs

Any additional keywords passed into this function will be passed on to the astropy test runner. This allows use of test-related functionality implemented in later versions of astropy without explicitly updating the package template.

unwrap_phase

shampoo.unwrap_phase(*reconstructed_wave*, *seed*=42)

2D phase unwrap a complex reconstructed wave.

Essentially a wrapper around the unwrap_phase function. The output will be type float64.

Parameters

reconstructed_wave : ndarray

Complex reconstructed wave

seed : float (optional)

Random seed, optional.

Returns

ndarray

Unwrapped phase image

4.1.2 Classes

<code>FFT(shape, float_precision, complex_precision)</code>	Convenience wrapper around <code>pyfftw.builders.fft2</code> .
<code>Hologram(hologram[, crop_fraction, ...])</code>	Container for holograms and methods to reconstruct them.
<code>ReconstructedWave(reconstructed_wave)</code>	Container for reconstructed waves and their intensity and phase arrays.

FFT

class shampoo.FFT(*shape*, *float_precision*, *complex_precision*, *threads*=2)

Bases: `object`

Convenience wrapper around `pyfftw.builders.fft2`.

Parameters

shape : tuple

Shape of the arrays which you will take the Fourier transforms of.

float_precision : dtype

complex_precision : dtype

threads : int, optional

This FFT implementation uses multithreading, with two threads by default.

Methods Summary

<code>fft2(array)</code>	2D Fourier transform.
<code>ifft2(array)</code>	Inverse 2D Fourier transform.

Methods Documentation

fft2(array)

2D Fourier transform.

Parameters

array : `ndarray` (real)

Input array

Returns

ft_array : `ndarray` (complex)

Fourier transform of the input array

ifft2(array)

Inverse 2D Fourier transform.

Parameters

array : `ndarray`

Input array

Returns

ift_array : `ndarray`

Inverse Fourier transform of input array

Hologram

class shampoo.**Hologram**(*hologram*, *crop_fraction*=None, *wavelength*=4.05e-07, *rebin_factor*=1, *dx*=3.45e-06, *dy*=3.45e-06, *threads*=2)

Bases: `object`

Container for holograms and methods to reconstruct them.

Parameters

hologram : `ndarray`

Input hologram

crop_fraction : float

Fraction of the image to crop for analysis

wavelength : float [meters]

Wavelength of laser

rebin_factor : int

Rebin the image by factor `rebin_factor`. Must be an even integer.

dx : float [meters]

Pixel width in x-direction (unbinned)

dy : float [meters]

Pixel width in y-direction (unbinned)

Notes

Non-square holograms will be cropped to a square with the dimensions of the smallest dimension.

Methods Summary

<code>apodize(array[, alpha])</code>	Force the magnitude of an array to go to zero at the boundaries.
<code>Hologram.detect_specimens</code>	
<code>fourier_peak_centroid(fourier_arr[, ...])</code>	Calculate the centroid of the signal spike in Fourier space near the frequencies of the real image.
<code>Hologram.fourier_trans_of_impulse_resp_func</code>	
<code>from_tif(hologram_path, **kwargs)</code>	Load a hologram from a TIF file.
<code>get_digital_phase_mask(psi[, plots])</code>	Calculate the digital phase mask (i.e.
<code>real_image_mask(center_x, center_y, radius)</code>	Calculate the Fourier-space mask to isolate the real image
<code>reconstruct(propagation_distance[, ...])</code>	Reconstruct the wave at <code>propagation_distance</code> .
<code>Hologram.reconstruct_multithread</code>	

Methods Documentation

apodize(*array*, *alpha*=0.075)

Force the magnitude of an array to go to zero at the boundaries.

Parameters

array : `ndarray`

Array to apodize

alpha : float between zero and one

Alpha parameter for the Tukey window function. For best results, keep between 0.075 and 0.2.

Returns

apodized_arr : `ndarray`

Apodized array

fourier_peak_centroid(*fourier_arr*, *mask_radius*=None, *margin_factor*=0.1, *plot*=False)

Calculate the centroid of the signal spike in Fourier space near the frequencies of the real image.

Parameters

fourier_arr : `ndarray`

Fourier-transform of the hologram

margin_factor : int

Fraction of the length of the Fourier-transform of the hologram to ignore near the edges, where spurious peaks occur there.

plot : bool

Plot the peak-centroiding visualization of the fourier transform of the hologram? Default is False.

Returns

pixel : `ndarray`

Pixel at the centroid of the spike in Fourier transform of the hologram near the real image.

classmethod `from_tif(hologram_path, **kwargs)`

Load a hologram from a TIF file.

This class method takes the path to the TIF file as the first argument. All other arguments are the same as `Hologram`.

Parameters

hologram_path : str

Path to the hologram to load

get_digital_phase_mask(psi, plots=False)

Calculate the digital phase mask (i.e. reference wave), as in Colomb et al. 2006, Eqn. 26 [R3].

Fit for a second order polynomial, numerical parametric lens with least squares to remove tilt, spherical aberration.

Parameters

psi : `ndarray`

The product of the Fourier transform of the hologram and the Fourier transform of impulse response function

plots : bool

Display plots after calculation if `True`

Returns

phase_mask : `ndarray`

Digital phase mask, used for correcting phase aberrations.

real_image_mask(center_x, center_y, radius)

Calculate the Fourier-space mask to isolate the real image

Parameters

center_x : int

x centroid [pixels] of real image in Fourier space

center_y : int

y centroid [pixels] of real image in Fourier space

radius : float

Radial width of mask [pixels] to apply to the real image in Fourier space

Returns

mask : `ndarray`

Binary-valued mask centered on the real-image peak in the Fourier transform of the hologram.

reconstruct(propagation_distance, plot_aberration_correction=False, plot_fourier_peak=False, cache=False)

Reconstruct the wave at `propagation_distance`.

If `cache` is `True`, the reconstructed wave will be cached onto the `Hologram` object for quick retrieval.

Parameters

propagation_distance : float

Propagation distance [m]

plot_aberration_correction : bool

Plot the aberration correction visualization? Default is False.

plot_fourier_peak : bool

Plot the peak-centroiding visualization of the fourier transform of the hologram? Default is False.

cache : bool

Cache reconstructions onto the hologram object? Default is False.

Returns

reconstructed_wave : [ReconstructedWave](#)

The reconstructed wave.

ReconstructedWave

class shampoo.**ReconstructedWave**(*reconstructed_wave*)

Bases: [object](#)

Container for reconstructed waves and their intensity and phase arrays.

Attributes Summary

intensity	Reconstructed intensity (ndarray , real)
phase	Reconstructed, unwrapped phase (ndarray , real)
reconstructed_wave	Reconstructed wave (ndarray , complex)

Methods Summary

plot ([phase , intensity , all , cmap])	Plot the reconstructed phase and/or intensity images.
---	---

Attributes Documentation

intensity

Reconstructed intensity ([ndarray](#), real)

phase

Reconstructed, unwrapped phase ([ndarray](#), real)

Phase unwrapping comes from `unwrap_phase`.

reconstructed_wave

Reconstructed wave ([ndarray](#), complex)

Methods Documentation

plot(*phase=False, intensity=False, all=False, cmap=<matplotlib.colors.LinearSegmentedColormap object>*)

Plot the reconstructed phase and/or intensity images.

Parameters**phase** : bool

Toggle unwrapped phase plot. Default is False.

intensity : bool

Toggle intensity plot. Default is False.

all : bool

Toggle unwrapped phase plot and . Default is False.

cmap : [Colormap](#)

Matplotlib colormap for phase and intensity plots.

Returns


fig : [Figure](#)

Matplotlib figure object

ax : [Axes](#)

Matplotlib axis object

4.1.3 Class Inheritance Diagram



Part III

Authors

- Brett Morris (UW)

Bibliography

[R4] <https://www.osapublishing.org/oe/abstract.cfm?uri=oe-14-13-5895>

[R3] <http://www.ncbi.nlm.nih.gov/pubmed/16512526>

S

shampoo, [17](#)

A

apodize() (shampoo.Hologram method), 24

C

cluster_focus_peaks() (in module shampoo), 17

create_hdf5_archive() (in module shampoo), 18

F

FFT (class in shampoo), 22

fft2() (shampoo.FFT method), 23

fftshift() (in module shampoo), 18

find_focus_plane() (in module shampoo), 18

fourier_peak_centroid() (shampoo.Hologram method), 24

from_tif() (shampoo.Hologram class method), 25

G

get_digital_phase_mask() (shampoo.Hologram method),
25

glue_focus() (in module shampoo), 19

H

Hologram (class in shampoo), 23

I

ifft2() (shampoo.FFT method), 23

intensity (shampoo.ReconstructedWave attribute), 26

L

locate_specimens() (in module shampoo), 19

O

open_hdf5_archive() (in module shampoo), 20

P

phase (shampoo.ReconstructedWave attribute), 26

plot() (shampoo.ReconstructedWave method), 26

R

real_image_mask() (shampoo.Hologram method), 25

reconstruct() (shampoo.Hologram method), 25

reconstructed_wave (shampoo.ReconstructedWave
attribute), 26

ReconstructedWave (class in shampoo), 26

S

save_scaled_image() (in module shampoo), 20

shampoo (module), 17

T

test() (in module shampoo), 21

U

unwrap_phase() (in module shampoo), 22